INSTRUMENTS USED IN OPEN CASTRATION

Photo of Instrument **Key Points for Use** Video Name **Emasculators** Simultaneously cut and crush (Size ranges from spermatic cord 9",12",14", 16" and 19") Crushing part of tool should be placed close to the body Nut of tool should be in direction of "nut" (testes) of animal Reimer Emasculator -The two crushing surfaces 12.5" allows for better haemostasis (2 crushing and 1 cutting Cutting is performed at a later section) stage using separate handle of the device

Serra Emasculator	 Instrument draws the spermatic cord concentrically into the jaws, preventing slippage Haemostasis is achieved by compression, stretching and tearing of tissues Crushing and transection of the spermatic cord is accomplished simultaneously through one single closing movement of the jaws
Newberry Knife	 Force testicles up and insert blade through scrotum and below testicles Close knife and pull downwards and slice the scrotum Remove testicles manually using emasculator/Henderson castration tool

Henderson Tool



- Attached to a cordless drill
- Twisting action closes the severed cord
- Minimizes blood loss and risk of haemorrhage through the severed cord

https://ww w.youtube. com/watch ?v=_w4onj F_Mu0