
## STANDARD HYDROGEN ELECTRODE


A standard hydrogen electrode (SHE) consists of a platinised platinum electrode immersed in a solution of 1 mol dm<sup>-3</sup> H<sup>+</sup> ions. Hydrogen gas at a pressure of 1 atm and 25 °C is bubbled over the platinum electrode.



The platinum electrode has a large surface area so that it can establish an equilibrium between  $H_2(g)$  and  $H^+$  (aq) as rapidly as possible.


- Standard reduction potential (or standard electrode potential) is the voltage associated with a reduction reaction at an electrode when all solutions are 1 M and all gases are at 1 atm.
- It is impossible to obtain the standard electrode potential of a single half cell because only differences in potentials are measureable.
- Standard reduction potential (S.R.P) values or E°<sub>red</sub> of half-cells are referenced to a Standard Hydrogen Electrode (SHE).
- The standard reduction potential for hydrogen is 0 V:
   2H⁺ (aq, 1M) + 2 e⁻ → H₂ (g, 1 atm)
- SHE acts as either cathode or anode when connected to other half-cells.

## Determination of Eo<sub>red</sub> of Zn<sup>2+</sup>(aq)/Zn(s) half-cell



Zn<sup>2+</sup>(aq)/Zn(s) half-cell is is anode and SHE is cathode.

## Determination of Eo<sub>red</sub> of Cu<sup>2+</sup>(aq)/Cu(s) half-cell



SHE is anode and Cu<sup>2+</sup>(aq)/Cu(s) half-cell is cathode.

E°<sub>red</sub> of Cu²+(aq)/Cu half-cell and Zn²+(aq)/Zn half-cell can be calculated by substituting the cell voltage (or standard cell potential) measured by the voltmeter into the following formula:

$$E_{cell}^{o} = E_{cathode}^{o} - E_{anode}^{o}$$

$$+ 0.76V = E_{H^{+}/H_{2}}^{o} - E_{Zn^{2+}/Zn}^{o}$$

$$+ 0.76V = 0.00 \text{ V} - E_{Zn^{2+}/Zn}^{o}$$

$$+ E_{Zn^{2+}/Zn}^{o} = -0.76V$$

$$E_{cell}^{o} = E_{cathode}^{o} - E_{anode}^{o}$$

$$+ 0.34 \text{ V} = E_{Cu^{2+}/Cu}^{o} - E_{H^{+}/H_{2}}^{o}$$

$$+ 0.34 \text{ V} = E_{Cu^{2+}/Cu}^{o} - 0.00 \text{ V}$$

$$+ E_{Cu^{2+}/Cu}^{o} = + 0.34 \text{ V}$$

$$E_{H^+/H_2}^o = 0.00 \text{ V}$$

$$E_{Zn^{2+}/Zn}^{o} = -0.76V$$

$$E_{Cu^{2+}/Cu}^{o} = + 0.34 \text{ V}$$

Standard reduction potential (S.R.P) or standard electrode potential

$$2H^{+}(aq) + 2e^{-} \longrightarrow H_{2}(g) E^{\circ} = 0.00 V$$

$$Zn^{2+}(aq) + 2e^{-} \longrightarrow Zn(s) \quad E^{\circ} = -0.76 \text{ V}$$

$$Cu^{2+}(aq) + 2e^{-} \longrightarrow Cu(s) E^{\circ} = + 0.34 V$$